

## Ultra High-Speed Analysis of Alkylphenones

### Introduction

Alkylphenones are often used for evaluation of HPLC systems because of their reasonable retentions on ODS columns and easy separation with water/acetonitrile mobile phase.


This report outlines the ultra-high speed analysis of Alkylphenones analysis by Ultra High-performance Liquid Chromatography (UHPLC) using with a PDA detector that enables high-speed data acquisition of 100 spectra/sec.

### Experimental

Equipment	Conditions
Pump: X-LC 3180PU x2	Column: ZORBAX SB-C18 (3.0 mmID x 30 mmL, 1.8 $\mu$ m)
Degasser: X-LC 3080DG	Eluent A: Water
Mixer: X-LC 3180MX	Eluent B: Acetonitrile
Column oven: X-LC 3161CO	Gradient condition: (A/B), 0 min(40/60) $\rightarrow$ 0.15 min(5/95) $\rightarrow$ 0.30 min(5/95) $\rightarrow$ 0.35 min(40/60) 1 cycle; 1 min
Autosampler: X-LC 3159AS	Flow rate: 3.8 mL/min
Detector: X-LC 3110MD	Column temp.: 40°C
	Wavelength: 200-300 nm
	Injection volume: 5 $\mu$ L
	Standard sample: Acetanilide, Acetophenone, Propiophenone, Butyrophenone, Benzophenone, Hexanophenone, Valerophenone, Heptanophenone, Octanophenone (50 $\mu$ g/mL each)

### Results

Fig. 1 shows the 3D Chromatogram of the alkylphenones standard mixture containing 9 components. Alkylphenones were successfully separated and detected within 12 seconds.



**Fig. 1.** 3D Chromatogram of standard alkylphenones. 1: Acetanilide, 2: Acetophenone, 3: Propiophenone, 4: Butyrophenone, 5: Benzophenone, 6: Valerophenone, 7: Hexanophenone, 8: Heptanophenone, 9: Octanophenone.