

Extraction of Methylmalonic Acid from Serum Using ISOLUTE® SAX Prior to LC-MS/MS Analysis




Figure 1. Structures of methylmalonic acid (MMA) and succinic acid (SA).

Introduction

Methylmalonic acid (MMA) in serum is measured to help diagnose a number of disorders, primarily Vitamin B₁₂ deficiency. This application note describes a simple, effective protocol for the extraction of methylmalonic acid (MMA) from serum using ISOLUTE® SAX strong anion exchange solid phase extraction plates, demonstrating high, reproducible analyte recoveries with low protein and phospholipid content in the extracts. The well known isobaric interference, succinic acid, is chromatographically separated to allow accurate quantitation of the MMA.

ISOLUTE SAX is a silica-based strong anion exchange sorbent used to extract acidic analytes, such as MMA, from aqueous samples. Very clean extracts are possible due to the purely anion exchange retention mechanism available using this silica based sorbent. This minimizes the co-extraction of endogenous interferences through hydrophobic (or non-polar) retention.

Analytes

MMA and MMA-¹³C₄ as internal standard.

Sample Preparation Procedure

Format:

ISOLUTE® SAX 25 mg fixed well plate, part number 500-0025-P01

Sample Pre-treatment

To serum (100 µL), add 10 µL of internal standard (10 ng/µL). Allow to stand for ~1 hour to allow binding to occur. Add HPLC grade water (190 µL) and vortex.

Conditioning

Condition each well with methanol (500 µL).

Equilibration

Equilibrate each well with HPLC grade water (500 µL).

Sample Loading

Load 300 µL of pre-treated sample.

Wash 1

Elute interferences with HPLC grade water (500 µL).

Wash 2

Elute interferences with methanol (500 µL).

Analyte Elution

Elute analytes with 2% formic acid in acetonitrile (600 µL) into a collection plate.

Post Extraction

Dry the extract in a stream of air or nitrogen using a SPE Dry (40 °C at 40 L/min) or TurboVap (40 °C at 1.0 bar).

Reconstitution

Add 100 µL of 0.4% formic acid (aq), seal with a plate mat and vortex for 30 seconds.

UPLC Conditions

Instrument

Waters ACQUITY I Class UPLC equipped with a flow through needle (15 µL)

Column

Gemini 3 µm C18 (100 x 3 mm id)

Mobile Phase

A: 0.4% formic acid (aq)

B: 0.4% formic acid in methanol

Flow Rate

0.6 mL/min

Table 1. Gradient Conditions - numerical.

Step	%A	%B	Curve
0	100	0	1
1	100	0	6
2.5	98	2	6
3	100	0	11

Curve 6: Linear Gradient

Curve 11: Conditions in line initiated immediately once time reached. i.e. 0% B resumed at 3 minutes.




Figure 2. Gradient Conditions - graphical

Injection Volume

10 µL

Sample Temperature

20 °C

Column Temperature

50 °C

MS Conditions

Instrument

Waters XEVO TQS triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis.

Desolvation Temperature:

500 °C

Ion Source Temperature:

150 °C

Negative ions were acquired in the multiple reaction monitoring (MRM) mode:

Table 2. MRM Conditions.

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
MMA	116.9 > 72.9	30	9
MMA- ¹³ C ₄	121.0 > 76.0	30	9

Results

Chromatography

Good separation was achieved between MMA and the isobaric interference succinic acid. **Figure 3.** shows a chromatogram of serum spiked with 10 ng/mL MMA and the baseline raised to 10:1 signal:noise, indicating an approximate lower limit of quantitation.




Figure 3. Chromatogram of ¹³C₄ MMA (top) at 100 ng/mL and MMA (bottom) at 10 ng/mL (~0.085 µMol/L) with x10 signal:noise indicator for the latter.

Recovery

Serum free of MMA was spiked at 250 ng/mL (~2.11 µMol/L). High reproducible recoveries >90% and corresponding RSDs of <10% were demonstrated. Typical recovery data is shown in Figure 4.




Figure 4. Chart demonstrating typical MMA recovery using the optimized ISOLUTE® SAX method.

Calibration Curves

Good linearity was observed over the range 10 – 2000 ng/mL (~0.085 - ~16.949 µMol/L). Figure 5 shows the coefficient of determination r^2 for the optimized method. In addition,

Compound name: MMA
Correlation coefficient: $r = 0.999977$, $r^2 = 0.999953$
Calibration curve: $0.597487 \times x + 4.97752$
Response type: Internal Std (Ref 2), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Include, Weighting: Null, Axis trans: None




Figure 5. Calibration line of spiked serum extracted with the optimized protocol.

Table 3. Calculated MMA concentrations.

Chromsystems Calibration Level	Set Value (ng/mL)	Calculated Value (ng/mL)
Calibrator 1	13.7	12.4
Calibrator 2	28.4	27.0
Calibrator 3	51.3	52.8

Additional Notes

Processing Guidelines

- » 96-well SPE plates were processed using a Biotage® PRESSURE+96 Positive Pressure Manifold at a pressure of 1–2 psi

Solvent Composition and Preparation Instructions

- » All solvents were HPLC grade.
- » 2% formic acid in acetonitrile: Add 200 µL concentrated formic acid to 9.8 mL of HPLC grade acetonitrile. Mix.
- » 0.4% formic acid (aq): Add 200 µL concentrated formic acid to 49.8 mL of HPLC grade water. Mix.
- » 0.4% formic acid in methanol: Add 200 µL concentrated formic acid to 49.8 mL of HPLC grade methanol. Mix.

Ordering Information

Part Number	Description	Quantity
500-0025-P01	ISOLUTE® SAX 25 mg Fixed Well Plate	1
121-5203	Collection plate, 2 mL, square	50
PPM-96	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
SD-9600-DHS-EU	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
SD-9600-DHS-NA	Biotage® SPE Dry Sample Concentrator System 100/120 V	1
C103264	TurboVap® 96, Evaporator 220/240V	1
C103263	TurboVap® 96, Evaporator 100/120V	1

Appendix

Biotage® Extrahera™ Settings

The method described in this application note was automated on the Biotage® Extrahera™, using ISOLUTE® SAX 25 mg SPE plates. Total time taken to process a full 96-well plate was 43 minutes. Method performance was comparable.


This appendix contains the software settings required to configure Extrahera to run this method. An importable electronic copy of this method for Extrahera can be downloaded from www.biotage.com

Biotage® Extrahera™ Data

Analyte	Methylmalonic Acid
Recovery (n=8) at 100 ng/mL	92.4%
%RSD	2.5
Linearity (r^2)	0.999*
LLOQ	<10 ng/ mL

*Note: Linearity experiments on Extrahera were run using 3PLUS1® Multilevel Plasma Calibrator Set Methylmalonic acid (Chromsystems Instruments and Chemicals GmbH). Manual processing using these standards gave linearity (r^2) of 0.999.

Data (manual processing) in the application note was generated using 'in house' spiked MMA free serum from Golden West Biologicals, Inc.



Method Name: MMA ISOLUTE® SAX
Sample Plate/Rack: 2 mL x 96 well 200 µL
Extraction Media: ISOLUTE SAX 96 Well Plate

< Cancel Edit SPE Method - MMA ISOLUTE SAX Save >

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Pretreatment: ☒ On Sample Pretreatment Conditioning Equilibration Load Wash (2) Elution

Conditioning: ☒ On Sample type: Aqueous sample Method comment:

Equilibration: ☒ On Starting sample volume in plate/rack (µL): 150

Load: ☒ On

Wash: ☒ On

Elution: ☒ On

Settings

"Sample" Tab

Sample Type:

Aqueous Sample

Starting Sample Volume (µL):

150

Method Comment:

Screenshot

Edit SPE Method - MMA ISOLUTE SAX

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Pre-treatment **Sample** **Pretreatment** **Conditioning** **Equilibration** **Load** **Wash (2)** **Elution**

Pre-treatment: ☒ On

Conditioning: ☒ On

Equilibration: ☒ On

Load: ☒ On

Wash: ☒ On

Elution: ☒ On

Number of steps: 1

Solvent: Water

Volume (µL): 300

Wait time (min): 0

Pause after last step? ☐ No

Dispose tips after each step? ☐ No

Settings

Pre-treatment	Activated
No. of steps	1
Pause after last step	No
Dispose tips after last step	No

Solvent
1 Water
2
3
4

	1	2	3	4
Volume (µL)	300			
Wait Time (min)	0			

Edit SPE Method - MMA ISOLUTE SAX

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Pre-treatment **Sample** **Pretreatment** **Conditioning** **Equilibration** **Load** **Wash (2)** **Elution**

Pre-treatment: ☒ On

Conditioning: ☒ On

Equilibration: ☒ On

Load: ☒ On

Wash: ☒ On

Elution: ☒ On

Number of steps: 1 Pressure (bar): 1.0

Solvent: Methanol

Volume (µL): 500 Collect in position: D (Waste)

Positive pressure time (s): 40 Advanced..

Repeat (number of times): 1 Pause after this step? ☐ No

Dispose tips after each step? ☐ No

Conditioning	Activated
No. of steps	1
Pressure (Bar)	1.0
Dispose tips after this step	No

Solvent
1 Methanol
2
3
4

	1	2	3	4
Volume (µL)	500			
Position	D			
Positive pressure time (s)	40			
Repeat	1			
Pause after this step	No			

Advanced Settings

< Cancel Edit SPE Method - MMA ISOLUTE SAX Save >

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Pre-treatment: ☒ On

Conditioning: ☒ On

Equilibration: ☒ On

Load: ☒ On

Wash: ☒ On

Elution: ☒ On

Number of steps: 1 Pressure (bar): 1.0

Solvent: Water

Volume (µL): 500 Collect in position: D (Waste)

Positive pressure time (s): 40 Advanced..

Repeat (number of times): 1 Pause after this step? ☐ No

Dispose tips after each step? ☐ No

Equilibration	Not Activated
No. of steps	1
Pressure (Bar)	1.0
Dispose tips after this step	No

Solvent	
1	Water
2	
3	
4	

	1	2	3	4
Volume (µL)	500			
Position	D			
Positive Pressure time (s)	40			
Repeat	1			
Pause after this step	No			

Advanced Settings

< Cancel Edit SPE Method - MMA ISOLUTE SAX Save >

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Pre-treatment: ☒ On

Conditioning: ☒ On

Equilibration: ☒ On

Load: ☒ On

Wash: ☒ On

Elution: ☒ On

Pressure (bar): 1.0 Pause after each load? ☐ No

Volume (µL): 300 Collect in position: D (Waste)

Positive pressure time (s): 80 Advanced..

Premix? ☒ Yes Number of times: 2

Load	
Pressure (Bar)	1.0
Pause after each load	No
Volume (µL)	300
Collect in position	D
Positive pressure time (s)	80
Premix	Yes
Number of times	2

Advanced Settings

Edit SPE Method - MMA ISOLUTE SAX

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Wash (2)

Pretreatment: ☒ On
 Conditioning: ☒ On
 Equilibration: ☒ On
 Load: ☒ On
 Wash: ☒ On
 Elution: ☒ On

Number of steps: 2 Pressure (bar): 1.0 Plate dry after last wash? ☒ Yes Plate dry time (s): 300
 Dispose tips after each step? ☐ No

1
 Solvent: Water Volume (µL): 500 Collect in position: D (Waste)
 Positive pressure time (s): 60 Repeat (number of times): 1 Pause after this step? ☐ No

2
 Solvent: Methanol Volume (µL): 500 Collect in position: D (Waste)
 Positive pressure time (s): 60 Repeat (number of times): 1 Pause after this step? ☐ No

Wash	Activated
No. of steps	2
Pressure (Bar)	1.0
Plate dry after last wash	Yes
Plate dry time (s)	300
Dispose tips after last step	No

Solvent
1 Water
2 Methanol
3
4

	1	2	3	4
Volume (µL)	500	500		
Position	D	D		
Positive pressure time (s)	60	60		
Repeat	1	1		
Pause after this step	No	No		

Advanced Settings

Edit SPE Method - MMA ISOLUTE SAX

Method name: MMA ISOLUTE SAX Sample plate/rack: 2mL x 96 well 200µL Extraction media: ISOLUTE SAX 96 Well Plate

Elution

Pretreatment: ☒ On
 Conditioning: ☒ On
 Equilibration: ☒ On
 Load: ☒ On
 Wash: ☒ On
 Elution: ☒ On

Number of steps: 1 Pressure (bar): 1.0 Plate dry after last elution? ☒ Yes Plate dry time (s): 30
 Dispose tips after each step? ☐ No

1
 Solvent: 2% Formic in MeCN Volume (µL): 600 Collect in position: A
 Positive pressure time (s): 80 Repeat (number of times): 1 Pause after this step? ☐ No

Elution	Activated
No. of steps	1
Pressure (Bar)	1.0
Plate dry after last elution	Yes
Plate dry time (s)	30
Dispose tips after each step	No

Solvent
1 2% Formic in MeCN
2
3
4

	1	2	3	4
Volume (µL)	600			
Position	A			
Positive pressure time (s)	80			
Repeat	1			
Pause after this step	No			

Advanced Settings

Solvent Properties

Solvent Description	
1	Water
2	Methanol
3	2% Formic in MeCN
4	
5	
6	
7	
8	
9	
10	



Solvent	1	2	3	4	5	6	7	8	9	10
Reservoir Type	Refillable			Non Refillable						
Capacity	N/A	N/A	N/A							
Aspiration flow rate (mL/min)	10	10	10							
Dispense flow rate (mL/min)	20	20	20							
Lower air gap flow rate (mL/min)	20	20	20							
Lower air gap volume (µL)	5	5	5							
Upper air gap flow rate (mL/min)	20	120	120							
Upper air gap volume (µL)	100	100	100							
Upper air gap dispense pause	300	300	300							
Conditioning?	Yes	Yes	Yes							
Conditioning number of times	2	3	3							
Conditioning flow rate (mL/min)	20	20	20							
Chlorinated	No	No	No							
Serial dispense	No	No	No							

< Cancel		Edit Sample - Aqueous sample	Save >
Sample Sample name Aqueous sample		Air Gap Lower air gap flow rate (mL/min) 20	
Sample description Default settings for aqueous		Lower air gap volume (µL) 5	
Aspiration flow rate (mL/min) 10		Upper air gap flow rate (mL/min) 120	
Dispense flow rate (mL/min) 20		Upper air gap volume (µL) 100	
		Upper air gap dispense pause (ms) 300	

"Sample" Screen

Sample name	Aqueous sample
Sample description	Aqueous sample
Aspiration flow rate (mL/min)	10
Dispense flow rate (mL/min)	20
Lower air gap flow rate (mL/min)	20
Lower air gap volume (µL)	5
Upper air gap flow rate (mL/min)	120
Upper air gap volume (µL)	100
Upper air gap dispense pause	300

< Cancel		Edit Extraction Media - ISOLUTE SAX 96 Well PL...	Save >
Extraction Media Name ISOLUTE SAX 96 Well Plate		Pipetting Height Solvent dispensation height (mm) -125.0	
Manufacturer Biotage		Sample dispensation height (mm) -135.0	
Part number 500-0025-P01		Aspiration height (mm) -135.0	
Sorbent load (mg) 25		Tune Pipetting Heights...	
Capacity volume (µL) 1000			
Format 96			
Comment 			

"Extraction Media" Screen

Name	ISOLUTE SAX
Manufacturer	Biotage
Part number	500-0025-P01
Sorbent load (mg)	25
Capacity volume (µL)	1000
Format	96
Comment	
Solvent dispensation height (mm)	-125.0
Sample dispensation height (mm)	-135.0
Aspiration height (mm)	-135.0

< Cancel		Edit Sample Plate/Rack - 2mL x 96 well 200µL	Save >
Sample Plate/Rack Name 2mL x 96 well 200µL		Pipetting Height Aspiration height (mm) -161.0	
Capacity volume (µL) 1800		Pretreatment dispensation height (mm) -153.0	
Format 96		Tune Pipetting Heights...	

"Sample Plate/Rack" Screen

Name	2 mL Sample x 96 well 200 µL
Capacity volume (µL)	1800
Format	96
Aspiration height (mm)	-161.0
Pre-treatment dispensation height (mm)	-153.0

< Cancel	Edit Pipette Tip - 1000 µL Biotage tip	Save >
Pipette Tip Name 1000 µL Biotage tip Manufacturer Biotage Part number 414141 Capacity (µL) 1000 Length (mm) 95		

"Pipette tip" Screen

Name	1000 µL Biotage Tip
Manufacturer	Biotage
Part number	414141
Capacity (µL)	1000
Length (mm)	95

Additional Information

In this automated method, 150 µL of pre-spiked (IS) serum sample is mixed with 300 µL of water during the pre-treatment step. This gives a total volume of 450 µL, from which 300 µL is loaded.

EUROPE

Main Office: +46 18 565900
Toll Free: +800 18 565710
Fax: +46 18 591922
Order Tel: +46 18 565710
Order Fax: +46 18 565705
order@biotage.com
Support Tel: +46 18 56 59 11
Support Fax: + 46 18 56 57 11
eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
Support Tel: +1 800 446 4752
Outside US: +1 704 654 4900
us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
Fax: +86 21 2898 6153
cn_order@biotage.com
cn-1-pointsupport@biotage.com

To locate a distributor,
please visit our website
www.biotage.com

Part Number: AN849.V.1

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.