


# Extraction of Methylmalonic Acid from Serum Using ISOLUTE® SLE+ Prior to LC-MS/MS Analysis



**Figure 1.** Structures of methylmalonic acid (MMA) and succinic acid (SA).

## Introduction

Methylmalonic acid (MMA) in serum is measured to help diagnose a number of disorders, primarily Vitamin B12 deficiency. This application note describes the extraction of methylmalonic acid from serum, prior to LC-MS/MS analysis. An effective and efficient ISOLUTE® SLE+ protocol optimized for the 200 µL capacity 96-well plate format is used. The simple sample preparation procedure delivers clean extracts and analyte recoveries greater than 80% with RSDs below 10%.

ISOLUTE® SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation.

## Analytes

MMA and MMA-<sup>13</sup>C<sub>4</sub> as internal standard.

## Sample Preparation Procedure

### Format:

ISOLUTE® SLE+ 200 µL Supported Liquid Extraction plate, part number 820-0200-P01

### Sample Pre-treatment

To serum (100 µL), add 10 µL of ISTD (10 ng/µL). Allow to equilibrate and add 4.6M formic acid (aq) (100 µL). Mix.

### Sample Loading

Load the pre-treated serum (200 µL) onto the plate and apply a pulse of vacuum or positive pressure (3–5 seconds) to initiate flow. Allow the sample to absorb for 5 minutes.

### Analyte Extraction

Apply MTBE (750 µL) and allow to flow under gravity for 5 minutes. Apply vacuum or positive pressure (5–10 seconds) to pull through any remaining extraction solvent.

### Post Elution

Evaporate the extract in a stream of air or nitrogen using a SPE Dry (40 °C, 20 to 40 L/min). Addition of 2 µL of ethylene glycol prior to evaporation can help reduce analyte losses due to volatility (optional).

### Reconstitution

Reconstitute extracts with 100 µL of 0.4% formic acid (aq).

## UPLC Conditions

### Instrument

Waters ACQUITY I Class UPLC equipped with a flow through needle (15 µL)

### Column

Gemini 3 µm C18 (100 x 3 mm id)

### Mobile Phase

A: 0.4% formic acid (aq)

B: 0.4% formic acid in methanol

### Flow Rate


0.6 mL/min

**Table 1.** Gradient Conditions - numerical.

Step	%A	%B	Curve
0	100	0	1
1	100	0	6
2.5	98	2	6
3	100	0	11

**Curve 6:** Linear Gradient

**Curve 11:** Conditions in line initiated immediately once time reached. i.e. 100% resumed at 3 minutes.



**Figure 2.** Gradient Conditions - graphical

### Injection Volume

10 µL

### Sample Temperature

20 °C

### Column Temperature

50 °C

## MS Conditions

### Instrument

Waters XEVO TQS triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis.

### Desolvation Temperature:

500 °C

### Ion Source Temperature:

150 °C

Negative ions were acquired in the multiple reaction monitoring (MRM) mode:


**Table 2.** MRM Conditions.

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
MMA	116.9 > 72.9	30	9
MMA- <sup>13</sup> C <sub>4</sub>	121.0 > 76.0	30	9

## Results

### Chromatography


Good separation was achieved between MMA and the isobaric interference succinic acid. **Figure 3.** demonstrates a chromatogram of serum spiked with 10 ng/mL MMA and the baseline raised to 10:1 signal:noise, indicating an approximate lower limit of quantitation.



**Figure 3.** Chromatogram of <sup>13</sup>C<sub>4</sub> MMA (top) at 100 ng/mL and MMA (bottom) at 10 ng/mL (~0.085 µMol/L) with x10 signal:noise indicator for the latter.

## Recovery

Serum free of MMA was spiked at 250 ng/mL (~2.11 µMol/L). High reproducible recoveries >80% and corresponding RSDs of <10% were demonstrated. Typical recovery data is shown in **Figure 4**.




**Figure 4.** Chart demonstrating typical MMA recovery using the optimized ISOLUTE® SLE+ method.

## Calibration Curves

Good linearity was observed over the range 10 – 2000 ng/mL (~0.085 – ~16.949 µMol/L). **Figure 5** demonstrates the coefficient of determination  $r^2$  for the optimized method. In addition, commercial calibration samples from plasma matrix were extracted and their concentration values were evaluated against the in-house calibration line. Good agreement was reached and concentrations are summarized in **Table 3**.

Compound name: MMA  
 Correlation coefficient:  $r = 0.999823$ ,  $r^2 = 0.999667$   
 Calibration curve:  $0.588726 * x + 4.75283$   
 Response type: Internal Std (Ref 2), Area \* (IS Conc / IS Area)  
 Curve type: Linear, Origin: Include, Weighting: Null, Axis trans: None



**Figure 5.** Calibration line of spiked serum extracted with the optimized protocol.

**Table 3.** Calculated MMA concentrations.

Chromsystems Calibration Level	Set Value (ng/mL)	Calculated Value (ng/mL)
Calibrator 1	13.7	10.0
Calibrator 2	28.4	27.8
Calibrator 3	51.3	56.5

## Additional Notes

### Solvent composition and preparation instructions

- » All solvents were HPLC grade.
- » 4.6M (17.7%) formic acid (aq): Add 8.85 mL of concentrated formic acid to 41.15 mL of HPLC grade water.
- » 0.4% formic acid (aq): Add 200 µL concentrated formic acid to 49.8 mL of HPLC grade water.
- » 0.4% formic acid in methanol: Add 200 µL concentrated formic acid to 49.8 mL of HPLC grade methanol

## Ordering Information

Part Number	Description	Quantity
<b>820-0200-P01</b>	ISOLUTE® 200 µL Supported Liquid Extraction Plate	1
<b>121-5203</b>	Collection plate, 2 mL, square	50
<b>PPM-96</b>	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
<b>SD-9600-DHS-EU</b>	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
<b>SD-9600-DHS-NA</b>	Biotage® SPE Dry Sample Concentrator System 100/120 V	1
<b>C103264</b>	TurboVap® 96, Evaporator 220/240V	1
<b>C103263</b>	TurboVap® 96, Evaporator 100/120V	1

# Appendix

## Biotage® Extrahera™ Settings

The method described in this application note was automated on the Biotage® Extrahera™, using ISOLUTE® SLE+ Supported Liquid Extraction plates. Total time taken to process a full 96-well plate was 27 minutes. Method performance was comparable.


This appendix contains the software settings required to configure Extrahera to run this method. An importable electronic copy of this method for Extrahera can be downloaded from [www.biotage.com](http://www.biotage.com)

### Biotage® Extrahera™ Data

Analyte	Methylmalonic Acid
Recovery (n=8) at 100 ng/mL	74.8%
%RSD	3.2
Linearity (r <sup>2</sup> )	0.999*
LLOQ	<10 ng/ mL

\*Note: Linearity experiments on Extrahera were run using 3PLUS1® Multilevel Plasma Calibrator Set Methylmalonic acid (Chromsystems Instruments and Chemicals GmbH). Manual processing using these standards gave linearity (r<sup>2</sup>) of 0.999.

Data (manual processing) in the application note was generated using 'in house' spiked MMA free serum from Golden West Biologicals, Inc.



**Method Name:** MMA ISOLUTE® SLE+ 200  
**Sample Plate/Rack:** 2 mL sample plate, 96  
**Extraction Media:** ISOLUTE SLE+ 200 µL 96 Well Plate

< Cancel
Edit SLE Method - MMA ISOLUTE SLE+ 200
Save >

Method name  
MMA ISOLUTE SLE+ 200

Sample plate/rack  
2 mL sample plate, 96

Extraction media  
ISOLUTE SLE+ 200µL 96 well.

Pretreatment   
 Load   
 Elution

Sample  
 Sample type: aq sample  
 Starting sample volume in plate/rack (µL): 150

Method comment  
 SLE method for MMA in serum. Start with 0.15 mL of serum in a 2 mL collection block and run on an SLE plus 200 plate with MTBE as the elution solvent. The method uses a highly acidic pretreatment solvent, as a result it is recommended to use this in one of the non-refillable front reservoirs and that this is safely washed after use.

## Settings

### "Sample" Tab

**Sample Type:** Aqueous Sample

**Starting Sample Volume (µL):** 150

**Method Comment:**

SLE method for MMA in serum. Start with 0.15 mL of serum in a 2 mL collection block and run on an ISOLUTE SLE+ 200 plate with MTBE as the elution solvent. The method uses a highly acidic pretreatment solvent, as a result it is recommended to use this in one of the non-refillable front reservoirs and that this is safely washed after use.

## Screenshot

**Edit SLE Method - MMA ISOLUTE SLE+ 200**

Method name: MMA ISOLUTE SLE+ 200 | Sample plate/rack: 2 mL sample plate, 96 | Extraction media: ISOLUTE SLE+ 200µL 96 well

**Pre-treatment** | **Sample** | Pretreatment | Load | Elution

Pre-treatment:  On  
 Load:  On  
 Elution:  On

Number of steps: 1

Pause after last step?  No  No  
 Dispose tips after each step?  No  No

Solvent: Formic Acid 4.6 M  
 Volume (µL): 150  
 Wait time (min): 0

## Settings

Pre-treatment	Activated
No. of steps	1
Pause after last step	Yes
Dispose tips after last step	No

Solvent	
1	Formic Acid 4.6 M
2	
3	
4	

	1	2	3	4
Volume (µL)	150			
Wait Time (min)	0			

**Edit SLE Method - MMA ISOLUTE SLE+ 200**

Method name: MMA ISOLUTE SLE+ 200 | Sample plate/rack: 2 mL sample plate, 96 | Extraction media: ISOLUTE SLE+ 200µL 96 well

**Pre-treatment** | **Sample** | Pretreatment | **Load** | Elution

Pre-treatment:  On  
 Load:  On  
 Elution:  On

Volume (µL): 200 | Collect in position: D (Waste) | Pause after each load?  No

Air push time (s): 15  
 Wait time (min): 5

Premix?  Yes | Number of times: 5

Load	
Air Push Time(s)	15
Pause after each load	No
Volume (µL)	200
Collect in position	D
Air Push Time (s)	15
Wait time (min)	5
Premix	Yes
Number of times	5

**Edit SLE Method - MMA ISOLUTE SLE+ 200**

Method name: MMA ISOLUTE SLE+ 200 | Sample plate/rack: 2 mL sample plate, 96 | Extraction media: ISOLUTE SLE+ 200µL 96 wellbr.

**Elution**

Number of steps: 1

Air push after last elution?  Yes | Air push time (s): 45 | Dispose tips after each step?  No

Solvent: MTBE

Volume (µL): 750 | Collect in position: A

Wait time (min): 5 | Advanced...


Repeat (number of times): 1 | Pause after this step?  No

Elution	Activated
No. of steps	1
Pause after last step	No
Air push after last elution	Yes
Air push time (s)	45
Dispose tips after each step	No

Solvent
1 MTBE
2
3
4

	1	2	3	4
Volume (µL)	750			
Collect in position	A			
Wait time (min)	5			
Repeat	1			
Pause	No			

Advanced Settings
No



## Solvent Properties

Solvent Description	
1	Methyl Tert Butyl Ether
2	
3	
4	
5	
6	4.6 M formic acid in water
7	
8	
9	
10	



Solvent	1	2	3	4	5	6	7	8	9	10
<b>Reservoir Type</b>	<b>Refillable</b>					<b>Non Refillable</b>				
Capacity	N/A					N/A				
Aspiration flow rate (mL/min)	10					10				
Dispense flow rate (mL/min)	10					20				
Lower air gap flow rate (mL/min)	10					20				
Lower air gap volume (µL)	5					5				
Upper air gap flow rate (mL/min)	120					120				
Upper air gap volume (µL)	100					100				
Upper air gap dispense pause	300					300				
Conditioning?	Yes					Yes				
Conditioning number of times	2					2				
Conditioning flow rate (mL/min)	10					20				
Chlorinated	No					No				
Serial dispense	No					No				

< Cancel Edit Sample - Aqueous sample Save >

Sample	Air Gap
Sample name Aqueous sample	Lower air gap flow rate (mL/min) 20
Sample description Default settings for aqueous	Lower air gap volume (µL) 5
Aspiration flow rate (mL/min) 10	Upper air gap flow rate (mL/min) 120
Dispense flow rate (mL/min) 20	Upper air gap volume (µL) 100
	Upper air gap dispense pause (ms) 300

**"Sample" Screen**

Sample name	Aqueous sample
Sample description	Aqueous sample
Aspiration flow rate (mL/min)	10
Dispense flow rate (mL/min)	20
Lower air gap flow rate (mL/min)	20
Lower air gap volume (µL)	5
Upper air gap flow rate (mL/min)	120
Upper air gap volume (µL)	100
Upper air gap dispense pause	300

< Cancel Edit Extraction Media - ISOLUTE SLE+ 200µL 96... Save >

Extraction Media	Pipetting Height
Name ISOLUTE SLE+ 200µL 96 well p	Solvent dispensation height (mm) -125.0
Manufacturer Biotage	Sample dispensation height (mm) -130.5
Part number 820-0200-P01	Aspiration height (mm) -148.5
Sorbent load (mg) 0	<a href="#">Tune Pipetting Heights...</a>
Capacity volume (µL) 0	
Format 96	
Comment	

**"Extraction Media" Screen**

Name	ISOLUTE® SLE+ 200 µL 96 well plate
Manufacturer	Biotage
Part number	820-0200-P01
Sorbent load (mg)	0
Capacity volume (µL)	0
Format	96
Comment	N/A
Solvent dispensation height (mm)	-125.0
Sample dispensation height (mm)	-130.5
Aspiration height (mm)	-148.5

< Cancel Edit Sample Plate/Rack - 2 mL sample plate, 96 Save >

Sample Plate/Rack	Pipetting Height
Name 2 mL sample plate, 96	Aspiration height (mm) -162.0
Capacity volume (µL) 1800	Pretreatment dispensation height (mm) -128.0
Format 96	<a href="#">Tune Pipetting Heights...</a>

**"Sample Plate/Rack" Screen**

Name	2 mL Sample plate, 96
Capacity volume (µL)	1800
Format	96
Aspiration height (mm)	-162.0
Pre-treatment dispensation height (mm)	-128.0



< Cancel	Edit Pipette Tip - 1000 µL Biotage tip	Save >
<div style="border: 1px solid gray; padding: 5px; width: fit-content; margin: auto;"> <p><b>Pipette Tip</b></p> <p>Name 1000 µL Biotage tip</p> <p>Manufacturer Biotage</p> <p>Part number 414141</p> <p>Capacity (µL) 1000</p> <p>Length (mm) 95</p> </div>		

**"Pipette tip" Screen**

Name	1000 µL Biotage Tip
Manufacturer	Biotage
Part number	414141
Capacity (µL)	1000
Length (mm)	95

## Additional Information

In this automated method, 150 µL of pre-spiked (IS) serum sample is mixed with 150 µL of water during the pre-treatment step. This gives a total volume of 300 µL, from which 200 µL is loaded.

### EUROPE

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: + 46 18 56 57 11  
 eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 2898 6655  
 Fax: +86 21 2898 6153  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

To locate a distributor,  
 please visit our website  
 www.biotage.com

### Part Number: AN851.V.1

© 2016 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.