


# Extraction of 25-hydroxy Vitamin D from Serum Using ISOLUTE® PLD+ Prior to LC-MS/MS Analysis

This application note describes the extraction of 25-hydroxy vitamin D from serum, prior to LC-MS/MS analysis.



**Figure 1.** Structures of 25-hydroxy Vitamin D.

## Introduction

ISOLUTE® PLD+ Protein and Phospholipid Removal plates offer a substantial improvement in extract cleanliness compared to traditional protein precipitation techniques for bioanalytical sample preparation.

This application note describes a simple, effective ISOLUTE® PLD+ protocol for the extraction of 25-hydroxy vitamin D from serum, demonstrating high, reproducible analyte recoveries with low protein and phospholipid content in the extracts.

## Analytes

25-hydroxy vitamin D<sub>2</sub>, 25-hydroxy vitamin D<sub>3</sub> and d<sub>6</sub>-25-hydroxy vitamin D<sub>3</sub> as the internal standard.

## Sample Preparation Procedure

**Format:** ISOLUTE® PLD+ Protein and Phospholipid Removal plate, part number 918-0050-P01

### Sample Pre-treatment

Add 10 µL of ISTD (equivalent to 30 ng/mL) to the serum sample. Mix. Allow to stand for ~1 hour for binding to occur.

### Solvent Application

Apply 400 µL of Acetonitrile (MeCN) to each well of the ISOLUTE® PLD+ plate.

### Sample Application

Add 100 µL of serum with ISTD and mix thoroughly via repeat aspirate/dispense steps.

### Analyte Elution


Apply vacuum -0.2 bar or 3 psi positive pressure for approximately 5 minutes. For highly particulate laden samples increased pressure or vacuum conditions may be required.

### Post Extraction

Dry the extract in a stream of air or nitrogen using a Biotage® SPE Dry (40 °C at 40 L/min) or TurboVap® (40 °C at 1.0 bar).

### Reconstitution

Apply 100 µL of 30/70 2 mM Ammonium Formate, 0.1% formic acid aq/MeOH.



**Figure 2.** Typical ISOLUTE® PLD+ procedure

## UPLC Conditions

### Instrument

Waters Acquity UPLC (Waters Assoc., Milford, MA, USA)

### Column

ACE EXCEL 2 C18-PFP, 100 mm x 2.1 mm id 2 µm, (ACT, UK)

### Mobile Phase

A: 2 mM ammonium formate/0.1% formic acid (aq)

B: 2 mM ammonium formate/0.1% formic acid/MeOH

### Flow Rate

0.4 mL/min

**Table 1.** UPLC Gradient Conditions.

Time	%A	%B	Curve
0	25	75	1
3	0	100	6
4	25	75	11

**Curve 11:** Conditions in line initiated immediately once time passed. i.e. 25:75 resumed at 4 minutes.

**Curve 6:** Linear Gradient

### Injection Volume

15 µL (partial loop with overflow)

### Sample Temperature

20 °C

### Column Temperature

40 °C

**Note:** alternative UPLC conditions may be suitable. Check for good chromatographic separation of isobaric interferences to ensure accurate analyte quantitation

## Mass Spectrometry Conditions

### Instrument

Quattro Premier XE triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) equipped with an electrospray interface for mass analysis

### Desolvation Temperature

450 °C

### Ion Source Temperature

150 °C

### Collision Cell Pressure

3.76 e<sup>-3</sup> mbar

Positive ions acquired in the multiple reaction monitoring (MRM) mode:


**Table 2.** Quattro Premier XE MRM parameters (Qualifier ion details shown in parenthesis).

Analyte	MRM Transition	Cone V	Collision Energy eV
25-OH D <sub>2</sub>	395.5 > 269.5 (395.5 > 119.2)	30	18 26
25-OH D <sub>3</sub>	383.5 > 257.5 (383.5 > 107.2)	30	17 25
d6-25-OH D <sub>3</sub>	389.6 > 263.5	30	16

## Results

### Chromatography


Good chromatographic separation of 25-hydroxy vitamin D<sub>2</sub> and D<sub>3</sub> was achieved in less than 3 minutes, as shown in **Figure 3**.



**Figure 3.** Chromatographic separation of 25-hydroxy vitamin D<sub>2</sub> and D<sub>3</sub> from Chromsystems calibrated serum at 14.8 and 19.6 ng/mL respectively. ISTD at 30 ng/mL.

## Recovery

Serum and stripped serum was spiked at various concentrations from 2–100 ng/mL. High reproducible recoveries > 70% with corresponding RSDs < 10% were demonstrated. Typical recovery data is shown in **Figure 4**.




**Figure 4.** Recovery profile for 25-hydroxy vitamin D extracted at 50 ng/mL.

## PBS/BSA Calibration Curves


Calibration curves were generated using PBS/BSA spiked at concentrations from 1–100 ng/mL. Good coefficients of determination were obtained for 25-hydroxy vitamin D2 and D3 ( $r^2 > 0.99$ ).

Compound name: 25-OH-Vitamin D2 (1)  
Correlation coefficient:  $r = 0.999130$ ,  $r^2 = 0.998264$   
Calibration curve:  $0.657782 * x + -0.0205817$   
Response type: Internal Std (Ref 3), Area \* (IS Conc. / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



**Figure 5.** PBS/BSA Calibration Curve for 25-OH vitamin D2 constructed from 1–100 ng/mL.

Compound name: 25-OH-Vitamin D3 (1)  
Correlation coefficient:  $r = 0.998026$ ,  $r^2 = 0.996065$   
Calibration curve:  $1.38183 * x + 0.315245$   
Response type: Internal Std (Ref 3), Area \* (IS Conc. / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




**Figure 6.** PBS/BSA Calibration Curve for 25-OH vitamin D3 constructed from 1–100 ng/mL.

## Chromsystems Calibration Curves


Curves were also generated using calibrated serum standards (obtained from Chromsystems) spiked at concentrations from 0–69 ng/mL. Good coefficients of determination were obtained for 25-hydroxy vitamin D2 and D3 ( $r^2 > 0.99$ ).

Compound name: 25-OH-Vitamin D2 (1)  
Correlation coefficient:  $r = 0.999723$ ,  $r^2 = 0.999446$   
Calibration curve:  $0.459786 * x + 0.298189$   
Response type: Internal Std (Ref 3), Area \* (IS Conc. / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



**Figure 7.** Chromsystems Calibrated Serum Curve for 25-OH vitamin D2 constructed from 4–69 ng/mL.

Compound name: 25-OH-Vitamin D3 (1)  
Correlation coefficient:  $r = 0.998829$ ,  $r^2 = 0.99857$   
Calibration curve:  $1.58919 * x + -0.862662$   
Response type: Internal Std (Ref 3), Area \* (IS Conc. / IS Area)  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



**Figure 8.** Chromsystems Calibrated Serum Curve for 25-OH vitamin D3 constructed from 0–58 ng/mL.

## DEQAS External Quality Assessment Scheme

Final method testing was performed for 5 DEQAS serum samples extracted alongside the Chromsystems calibrators using the optimized method. The DEQAS criteria for acceptable performance is that at least 80% of results should fall within + or - 25% of the All Laboratory Trimmed Mean. Method performance is shown in **Table 3**. Units are quoted as ng/mL. All values fall within the accepted criteria.

**Table 3.** DEQAS 25-OH vitamin D results obtained using optimum method.

DEQAS Sample I.D.	DEQAS LC/MS Mean	ISOLUTE® PLD+
451	12.9	14.5
452	46.7	49.1
453	26.6	28.9
454	21.4	25.3
455	22.2	23.7

## Ordering Information

Part Number	Description	Quantity
<b>918-0050-P01</b>	ISOLUTE® PLD+ Fixed Well Plate	1
<b>121-9600</b>	Biotage® VacMater™-96 Sample Processing Manifold	1
<b>PPM-96</b>	Biotage® PRESSURE+ 96 Positive Pressure Manifold	1
<b>SD-9600-DHS-EU</b>	Biotage® SPE Dry Sample Concentrator System 220/240 V	1
<b>SD-9600-DHS-NA</b>	Biotage® SPE Dry Sample Concentrator System 100/120 V	1
<b>C103263</b>	TurboVap®96, Evaporator 100/120V	1
<b>C103264</b>	TurboVap® 96, Evaporator 220/240V	1

## Additional Notes

### Buffer Preparation

- 2 mM ammonium formate/0.1% formic acid (aq): Weigh 0.12612 g and dissolve in H<sub>2</sub>O. Add 1 mL of formic acid and make up to 1 L in H<sub>2</sub>O.
- 2 mM ammonium formate/0.1% formic acid/MeOH: Weigh 0.12612 g and dissolve in MeOH. Add 1 mL of formic acid and make up to 1 L in MeOH.

### Processing Conditions

Positive Pressure: Process at approximately 3 psi.

Vacuum Processing: Process at approximately -0.2 bar.

# Appendix

## Biotage® Extrahera™ Settings

The method described in this application note was automated on the Biotage® Extrahera™, using ISOLUTE PLD+ Protein and Phospholipid Removal plates. Method performance was comparable.

This appendix contains the software settings required to configure Extrahera to run this method.

An importable electronic copy of this method for Extrahera can be downloaded from [www.biotage.com](http://www.biotage.com)

### Biotage® Extrahera™ Data

Analyte	25-OH Vitamin D2	25-OH Vitamin D3
Recovery (n=8) at 50 ng/mL	92.0	81.1
%RSD	4.0	7.4
Linearity (r <sup>2</sup> )	0.999	0.999
LLOQ	<4 ng/mL	<4 ng/mL

**Method name:** AN842 Biotage® Extrahera™  
25OH Vit D PLD+

**Sample plate/rack:** 2 mL sample plate, 96

**Extraction Media:** ISOLUTE® PLD+ plate



## Screenshot

< Cancel    Edit PPT/PLD Method - 25OH Vit D PLD+    Save >

Method name: 25OH Vit D PLD+    Sample plate/rack: 2 mL sample plate, 96    Extraction media: PLD+ for 25 OH Vit D

Solvent Load:  On

Collection:  On

Sample type: PPT/PLD Sample

Starting sample volume in plate/rack (µL): 110

Method comment:

## Settings

### "Sample" Tab

**Sample Type:** PPT/PLD Sample

**Starting Sample Volume (µL):** 110

**Method Comment:**


Transfer 100 µL of plasma to a PLD plate containing 400 µL acetonitrile. A 100 µL volume is mixed 5 times before a 0.4 bar pressure is applied for 5 minutes

## Screenshot

## Settings

Solvent Load		Activated
Dispose tips		No
Solvent		
1	Acetonitrile	
		1
Volume (µL)	400	

Collection		Activated
Volume (µL)	100	
Premix	Yes	
Number of times	3	
Mix number of times	5	
Mix volume (µL)	100	
Wait time (min)	0	
Pressure (bar)	0.4	
Positive pressure time (s)	300	
Collect in position	A	



## Solvent Properties

Solvent Description	
1	Acetonitrile
2	
3	
4	
5	
6	
7	
8	
9	
10	



Solvent	1	2	3	4	5	6	7	8	9	10
Reservoir type	Refillable					Non Refillable				
Capacity	N/A									
Aspiration flow rate (mL/min)	10									
Dispense flow rate (mL/min)	20									
Lower air gap flow rate (mL/min)	20									
Lower air gap volume (µL)	5									
Upper air gap flow rate (mL/min)	120									
Upper air gap volume (µL)	100									
Upper air gap dispense pause	300									
Conditioning?	Yes									
Conditioning number of times	2									
Conditioning flow rate (mL/min)	20									
Chlorinated	No									
Serial dispense	No									

### Screenshot

**Edit Sample - PPT/PLD Sample**

<b>Sample</b> Sample name PPT/PLD Sample Sample description Sample for PPT/PLD method: Aspiration flow rate (mL/min) 50 Dispense flow rate (mL/min) 220	<b>Air Gap</b> Lower air gap flow rate (mL/min) 20 Lower air gap volume (µL) 5 Upper air gap flow rate (mL/min) 120 Upper air gap volume (µL) 700 Upper air gap dispense pause (ms) 1000
---	--

**"Sample" Screen**

Sample name	PLD/PPT Sample
Sample description	PLD/PPT Sample
Aspiration flow rate (mL/min)	50
Dispense flow rate (mL/min)	220
Lower air gap flow rate (mL/min)	20
Lower air gap volume (µL)	5
Upper air gap flow rate (mL/min)	120
Upper air gap volume (µL)	700
Upper air gap dispense pause	1000

**Edit Extraction Media - PLD+ for 25 OH Vit D**

<b>Extraction Media</b> Name PLD+ for 25 OH Vit D Manufacturer Biotage Part number 918-0050-P01 Sorbent load (mg) 0 Capacity volume (µL) 0 Format 96 Comment	<b>Pipetting Height</b> Solvent dispensation height (mm) -135.0 Sample dispensation height (mm) -135.0 Aspiration height (mm) -151.0 Tune Pipetting Heights...
---	---

**"Extraction Media" Screen**

Name	ISOLUTE® PLD+
Manufacturer	Biotage
Part number	918-0050-P01
Sorbent load (mg)	0
Capacity volume (µL)	0
Format	96
Comment	N/A
Solvent dispensation height (mm)	-135.0
Sample dispensation height (mm)	-135.0
Aspiration height (mm)	-151.0

**Edit Sample Plate/Rack - 2 mL sample plate, 96**

<b>Sample Plate/Rack</b> Name 2 mL sample plate, 96 Capacity volume (µL) 1800 Format 96	<b>Pipetting Height</b> Aspiration height (mm) -162.0 Pretreatment dispensation height (mm) -128.0 Tune Pipetting Heights...
---	---

**"Sample Plate/Rack" Screen**

Name	2 mL Sample Plate, 96
Capacity volume (µL)	1800
Format	96
Aspiration height (mm)	-162.0
Pre-treatment dispensation height (mm)	-128.0

< Cancel		Edit Pipette Tip - 1000 µL Biotage tip	Save >
<div style="border: 1px solid #ccc; padding: 5px; width: fit-content; margin: 0 auto;"> <p><b>Pipette Tip</b></p> <p>Name</p> <input type="text" value="1000 µL Biotage tip"/> <p>Manufacturer</p> <input type="text" value="Biotage"/> <p>Part number</p> <input type="text" value="414141"/> <p>Capacity (µL)</p> <input type="text" value="1000"/> <p>Length (mm)</p> <input type="text" value="95"/> </div>			

"Pipette tip" Screen	
<b>Name</b>	1000 µL Biotage Tip
<b>Manufacturer</b>	Biotage
<b>Part number</b>	414141
<b>Capacity (µL)</b>	1000
<b>Length (mm)</b>	95

**Additional Comments**

In this application the ISOLUTE® PLD+ extraction media default setting are edited such that the tip mixes the very top 100 µL of the 500 µL mixture. If the default ISOLUTE® PLD+ settings are used there is a greater risk of tip blockage during mixing.

**EUROPE**

Main Office: +46 18 565900  
 Toll Free: +800 18 565710  
 Fax: +46 18 591922  
 Order Tel: +46 18 565710  
 Order Fax: +46 18 565705  
 order@biotage.com  
 Support Tel: +46 18 56 59 11  
 Support Fax: + 46 18 56 57 11  
 eu-1-pointsupport@biotage.com

**NORTH & LATIN AMERICA**

Main Office: +1 704 654 4900  
 Toll Free: +1 800 446 4752  
 Fax: +1 704 654 4917  
 Order Tel: +1 704 654 4900  
 Order Fax: +1 434 296 8217  
 ordermailbox@biotage.com  
 Support Tel: +1 800 446 4752  
 Outside US: +1 704 654 4900  
 us-1-pointsupport@biotage.com

**JAPAN**

Tel: +81 3 5627 3123  
 Fax: +81 3 5627 3121  
 jp\_order@biotage.com  
 jp-1-pointsupport@biotage.com

**CHINA**

Tel: +86 21 2898 6655  
 Fax: +86 21 2898 6153  
 cn\_order@biotage.com  
 cn-1-pointsupport@biotage.com

To locate a distributor, please visit our website [www.biotage.com](http://www.biotage.com)

**Part Number: AN842.V.1**

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

