

ISOLUTE[®] SLE+ Supported Liquid Extraction Products

Simple Load-Wait-Elute Methodology

Supported Liquid Extraction

Achieving Simplicity and Success in Sample Preparation

ISOLUTE® SLE+

Supported Liquid Extraction Columns and Plates

ISOLUTE® SLE+ Supported Liquid Extraction products are designed to provide stress free extraction of analytes from biological fluids, using a simple **Load-Wait-Elute** methodology.


The supported liquid extraction (SLE) process is analogous to traditional liquid-liquid extraction (LLE) and utilizes the same water immiscible solvent systems for analyte extraction. Instead of shaking the two immiscible phases together, in SLE, the aqueous sample is immobilized on an inert support, and the organic phase flows through the support, eliminating problems such as emulsion formation and low analyte recoveries.

Methods with high analyte recoveries, and clean, protein and phospholipid free extracts, are easy to develop and automation is simple.

Select the correct ISOLUTE SLE+ product based on the volume of sample to be extracted (see **Table 1**). Extraction solvent volumes are also listed.

Table 1.
Recommended sample and elution volumes for ISOLUTE® SLE+ products.

| Product Description | Maximum Load Volume | Elution Protocol/Volume |
|-------------------------|---------------------|--------------------------|
| 200 µL 96-well plate | 200 µL | 1 x 1 mL |
| 400 µL 96-well plate | 400 µL | 2 x 900 µL or 3 x 700 µL |
| 1 mL 48-well plate | 1 mL | 5 x 1 mL |
| 200 µL Array well/plate | 200 µL | 2 x 600 µL |
| 400 µL Array well/plate | 400 µL | 3 x 750 µL |
| 400 µL column | 400 µL | 2 x 900 µL |
| 1 mL column | 1 mL | 2 x 2.5 mL |
| 2 mL column | 2 mL | 2 x 5 mL |
| 5 mL column | 5 mL | 3 x 8 mL |
| 10 mL column | 10 mL | 2 x 20 mL |



Typical ISOLUTE SLE+ procedure.

Supported Liquid Extraction

Achieving Simplicity and Success in Sample Preparation



96-Well Plates

| Part Number | Description | Qty. |
|--------------|---|------|
| 820-0200-P01 | ISOLUTE SLE+ 200 µL Supported Liquid Extraction Plate | 1 |
| 820-0400-P01 | ISOLUTE SLE+ 400 µL Supported Liquid Extraction Plate | 1 |

48-Well Plates

Extract 1 mL sample volumes in high throughput microplate format

| Part Number | Description | Qty. |
|--------------|---|------|
| 820-1000-Q01 | ISOLUTE SLE+ 1 mL Supported Liquid Extraction Plate (48-well) | 1 |

Columns

| Part Number | Description | Qty. |
|-------------|---|------|
| 820-0055-B | ISOLUTE SLE+ 400 µL Sample Volume | 50 |
| 820-0055-BG | ISOLUTE SLE+ 400 µL Sample Volume (Tablets) | 50 |
| 820-0140-C | ISOLUTE SLE+ 1 mL Sample Volume | 30 |
| 820-0140-CG | ISOLUTE SLE+ 1 mL Sample Volume (Tablets) | 30 |
| 820-0290-D | ISOLUTE SLE+ 2 mL Sample Volume | 20 |
| 820-0690-E | ISOLUTE SLE+ 5 mL Sample Volume | 20 |
| 820-1420-F | ISOLUTE SLE+ 10 mL Sample Volume | 16 |

Bulk Packs

To reduce packaging and improve workflow, ISOLUTE® SLE+ columns are available in bulk pack sizes

| Part Number | Description | Qty. |
|------------------|---|------|
| 820-0055-B-500 | ISOLUTE SLE+ 400 µL Sample Volume | 500 |
| 820-0055-BG-500 | ISOLUTE SLE+ 400 µL Sample Volume (Tablets) | 500 |
| 820-0140-C-1000 | ISOLUTE SLE+ 1 mL Sample Volume | 1000 |
| 820-0140-CG-1000 | ISOLUTE SLE+ 1 mL Sample Volume (Tablets) | 1000 |
| 820-0290-D-1000 | ISOLUTE SLE+ 2 mL Sample Volume | 1000 |

Array Wells


| Part Number | Description | Qty. |
|-------------|---------------------------------|------|
| 820-0200-T | ISOLUTE SLE+ 200 µL Array Wells | 100 |
| 820-0400-T | ISOLUTE SLE+ 400 µL Array Wells | 100 |

Accessories

| Part Number | Description | Qty. |
|-------------|---------------------------------|------|
| 121-5202 | Collection plate, 1 mL, square | 50 |
| 121-5203 | Collection plate, 2 mL, square | 50 |
| 121-5213 | Collection plate, 2 mL, round | 50 |
| 121-5210 | Collection plate, 5 mL, 48-well | 20 |

ISOLUTE® SLE+ User Guide

For further information download the ISOLUTE® SLE+ User Guide from www.biotage.com. Literature part number UI304.



Biotage are constantly developing new applications on ISOLUTE SLE+ products. Visit www.biotage.com for the latest information.