

# Extraction of Low Level Testosterone and Androstenedione From Human Serum Samples Using ISOLUTE® SLE+




Figure 1. Structure of Testosterone

## Introduction

This application note describes the extraction of testosterone and androstenedione from female patient serum samples using ISOLUTE SLE+ (96 well) plates.

This method has been optimized to extract testosterone and androstenedione from female clinical serum samples. The method has achieved very low limits of quantitation for testosterone and androstenedione 0.4 nmol/L and 0.9 nmol/L respectively. The recoveries for both analytes were greater than 90%.

ISOLUTE SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid liquid extraction formation, and significantly reduced sample preparation time.

## Analytes

Testosterone and androstenedione

## ISOLUTE SLE+ procedure

ISOLUTE SLE+ 400 Supported Liquid Extraction plate, part number 820-0400-P01.

**Sample pre treatment:** dilute human serum (200 µL) with 0.5 mol/L ammonium hydroxide (200 µL) then mix.

**Sample load:** load pre treated sample (400 µL) to plate followed by a pulse of vacuum to imitate flow and leave for 5 minutes.

**Analyte elution:** apply 500 µL of diethyl ether, wait five minutes to allow the solvent to soak, apply a short pulse of vacuum if solvent not fully absorbed. Apply a second 500 µL of diethyl ether; allow to soak for five minutes, apply a short pulse of vacuum if not fully absorbed. Apply a third 500 µL of Diethyl ether, allow to soak for 5 minutes and then apply another short pulse of vacuum.

**Post extraction:** evaporate the eluate to dryness and reconstitute with 400 µL of methanol:water (1:1, vol:vol).

**Additional information:** testosterone has an affinity to bind to plastic so the extracts were collected in glass tubes held in a 96 well collection plate.

Application Note

## HPLC conditions

**Instrument:** Acquity UPLC

**Column:** Acquity UPLC HSS C18 2.1 x 50 mm column, 1.8 µm particle (Waters)

**Mobile phase:** A= Water with 2 mmol/L ammonium acetate and 0.1% formic acid, B= Methanol with 2mmol/L ammonium acetate and 0.1% formic acid at a flow rate of 0.6 mL/min.

## Gradient

Time	%A	%B
0	50	50
0.20	50	50
1.50	30	70
2.25	2	98
2.75	30	70

**Injection volume:** 37.5 µL

## Mass spectrometry conditions

**Instrument:** Waters Premier XE


**Ion source temperature:** 120 °C

**Desolvation temperature:** 450 °C

**Table 1:** MRM transitions for Testosterone, Androstenedione, D<sub>2</sub> testosterone and D<sub>7</sub> androstenedione

Analyte	MRM transitions	Cone voltage	Collision energy
Testosterone	289.2>97.0	35	25
D2 Testosterone	291.1>98.9	35	23
Androstenedione	287.2>96.9	36	22
D7 Androstenedione	294>99.8	35	23

## Results:



**Figure 2.** Average analyte recoveries up to 100 nmol/L

**Table 2:** LLOQ of Testosterone and Androstenedione

Analyte	LLOQ nmol/L
Testosterone	0.4
Androstenedione	0.9

**References**

'Development and evaluation of a simple extraction method for the simultaneous measurement of female serum testosterone and androstenedione by UPLC-MSMS'  
M.Purcell, P.Reed. Department of Clinical Biochemistry, Salford Royal Hospital, M6 8HD, UK. Poster presented at IBMS September 2011.

**Ordering Information**

Part number	Description	Quantity
820-0400-P01	ISOLUTE SLE+ 400 µL Supported Liquid Extraction Plate	1

**NORTH AMERICA**

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
US-1-pointsupport@biotage.com

**EUROPE**

Main Office: +46 18 56 5900  
Fax: +46 18 59 1922  
Order Tel: +46 18 56 57 10  
Order Fax: +46 18 56 57 05  
order@biotage.com  
EU-1-pointsupport@biotage.com

**JAPAN**

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
JP-1-pointsupport@biotage.com

**CHINA**

Tel: +86 21 2898 6655  
Fax: +86 21 2898 6153  
CN\_order@biotage.com  
CN-1-pointsupport@biotage.com  
www.biotage.cn

**Distributors**

To locate a distributor please visit our  
Web site at [www.biotage.com](http://www.biotage.com)