


# Developing a Scalable Metals Removal Process



**Figure 1.** Metal scavenger (left). Target metal bound to metal scavenger (right).

This application note demonstrates the practical application of metal scavenging in the removal of palladium from a model Active Pharmaceutical Ingredient (API).

## Introduction


Biotage supported resin and silica based metal scavengers are effective reagents for the removal of trace Platinum Group Metals (PGMs) from catalyzed reactions to final APIs. They are applied in a variety of different industries, from pharmaceutical to fine chemical, from agrochemical to waste treatment.

In this study, metal scavenger kit K-MS-3 was used to first screen a contaminated compound and provide the basis for an Standard Operating Procedure (SOP) for larger scale metal scavenging. The kit provides metal scavengers in an easy to use screening evaluation format with example SOPs to enable rapid determination of a suitable metal scavenging candidate for an API mixture. The metal scavengers are industry proven and supported by a regulatory qualification support package.

A 1 g cartridge of metal scavenger is a good starting point to evaluate the product of a metal catalyzed reaction that used 1–2 mol % catalyst at 1 mmol scale. For larger/smaller reactions, this can be scaled accordingly.

## Results and Discussion

### Screening



**Scheme 1.** The initial scavenger screen.

The purpose of the initial screen was to quickly identify candidates that are not likely to be effective downstream. Based on our data below, carbon (our control experiment) and ISOLUTE® SCX-2 were least effective in this application, so Biotage® MP-TMT, ISOLUTE® Si-Thiol, ISOLUTE® Si-TMT and ISOLUTE® Si-Trisamine were chosen as potential scale-up candidates.

In our work, we were also looking for method robustness, so we eliminated ISOLUTE Si-Trisamine as a candidate as its basic character may limit the downstream scope of application (for example if an API contains acidic sites). Thus ISOLUTE Si-TMT became the primary screening candidate for our scale-up trials.



**Figure 2.** Convenient ready to use metal scavenging screening cartridges.

Scavenger 1 g/15 mL Fixed Bed	Time to Flow Through/s	Vol. Flow Rate mL/min	Approx. Linear Flow Rate cm/min	Pd Before ppm	Pd After ppm	Mass Recovery %	Pd Reduction %
Biotage® MP-TMT	333	1.8	1.4	500	10	99.3	98.0
ISOLUTE® Si-TMT	440	1.0	0.8	500	9	99.3	98.2
ISOLUTE® Si-Thiol	430	1.0	0.8	500	10	98.1	98.0
ISOLUTE® Si-Trisamine	371	1.4	1.0	500	9	99.8	98.2
ISOLUTE® SCX-2	418	1.4	1.0	500	469	91.6	6.2
Carbon	380	1.6	1.2	500	94	66.0	81.2

**Table 1.** Results of initial scavenger screen. Four possible candidates are identified, shortlisted to three.

ISOLUTE® Si-TMT	Vol. Flow Rate mL/min	Approx. Linear Flow Rate cm/min	Pd Before/ppm	Pd After/ppm	Pd Reduction %
<b>1 g/15 mL</b>	1.0	0.8	500	9	98.2
<b>50 g fixed bed</b>	15	0.8	25000	95	99.6

**Table 2.** Consistent fifty fold scalable results from screening to lab.

## Scale-Up

Based on initial screening results, a number of parameters were determined for further scale-up. The aspect ratio of the fixed bed was kept consistent with the 1 g/15 mL screening experiment by the use of a 50 g Biotage® SNAP cartridge, but maintaining aspect ratio is not always necessary for metal scavenging depending on the mechanism of action.

The data shows that the scavenging performance from screen to larger scale was reliable. We noted a slight improvement in scavenging on larger scale, which we attribute to our ability to more accurately make measurements on a larger scale. It may be that the small scale screening was more effective than we recorded. Based on these results; rationale to apply metal scavenging at even larger scale exists using linear flow rates corresponding to the 50 g flow style scavenging example.

### Scale Up Parameters

Metal Scavenger: ISOLUTE® Si-TMT, Mass of Fixed-Bed: 50 g, Dimensions: 80 mm x 40 mm (SNAP 50 g cartridge), Flow Mode: one pass, Equilibration Protocol: 3CV clean ethyl acetate, Volumetric Flow Rate: 15 mL/min, Linear Flow Rate: 0.8 cm/min.

## Mass Yields

In our screen, we compared Biotage metal scavengers with a classic carbon treatment. The carbon adsorption method removed some of the palladium but it was not the most effective. More significantly, our API analogue was lost to the carbon column (recovered yields were generally 66% compared to near quantitative from the metal scavengers) denoting agreement with consensus on the generic impact of carbon adsorption processing on mass yield and atom economy. In all cases, Biotage metal scavengers afforded higher yields of API.

## Equipment

Metal Scavengers, Biotage® MP-TMT, ISOLUTE® Si-TMT, ISOLUTE® Si-Thiol, ISOLUTE® SCx-2, ISOLUTE® Si-Trisamine. Format: 1 g/15 mL cartridge. Solvents: standard HPLC grade lab solvents. Other: Biotage® Gravity Rack p/n 123-2016 / 123-2019 (16 mm/ 19 mm dia tube rack).

## Detailed Methodology

Scavenger cartridges were equilibrated using 3–5 bed volumes of clean solvent (no need to collect eluent). 1 mL (corresponding to approx. 100 mg of dissolved API) was added to each cartridge on a Biotage® Gravity Rack and allowed to flow through under gravity at RT. [In some positions, the scavenger bed became discolored at the surface, which is usually an indication of metal scavenging]. The cartridges were each further washed with 2–4 mL of ethyl acetate. Since the scavenged metal is tightly bound, this wash step may be increased significantly to remove any slower running API from the cartridge. The eluent from each of the tubes was evaporated to dryness. The extract was weighed and analyzed for Pd content by ICP analysis.

## Ordering Information

Part Number	Description	Quantity
<b>121-0001</b>	PTFE Stopcock/Needle Unit	10
<b>121-0009</b>	Universal PTFE Stopcock	10
<b>K-MS-2</b>	Metal Scavenging Toolkit - Batch	1
<b>K-MS-3</b>	Metal Scavenging Screening Kit - SPE	1
<b>1201-0126-D</b>	Top Cap, 15 mL Colum	100
<b>123-2016</b>	Biotage® Gravity rack with 16 mm Collection Tube Rack	1
<b>123-2019</b>	Gravity Rack with 19 mm Collection Tube Rack	1

### EUROPE

Main Office: +46 18 565900  
Toll Free: +800 18 565710  
Fax: +46 18 591922  
Order Tel: +46 18 565710  
Order Fax: +46 18 565705  
order@biotage.com  
Support Tel: +46 18 56 59 11  
Support Fax: + 46 18 56 57 11  
eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
Support Tel: +1 800 446 4752  
Outside US: +1 704 654 4900  
us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 2898 6655  
Fax: +86 21 2898 6153  
cn\_order@biotage.com  
cn-1-pointsupport@biotage.com

### KOREA

Tel: + 82 31 706 8500  
Fax: + 82 31 706 8510  
korea\_info@biotage.com  
kr-1-pointsupport@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

### Part Number: AN878

© 2017 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.