



## APPLICATION NOTE

### PAH Analysis in Drinking Water using Disk SPE Automated with the SPE-DEX 5000

Alicia Cannon, Horizon Technology, Salem, NH USA

#### Key Words

PAH, Polycyclic aromatic hydrocarbons, USEPA 550.1

#### Introduction

Polycyclic aromatic hydrocarbons (PAHs) are ubiquitous pollutants introduced into the environment from anthropogenic activities such as combustion of fossil fuels, tobacco smoking and even during food preparation. These compounds strongly adsorb to sediment or particulate matter, causing them to slowly biodegrade while remaining present in drinking water sources; mostly as leachates from water storage and distribution systems. Exposure to PAHs over time can lead to increased health risks including developmental and reproductive issues. This class of compounds has also been identified as a carcinogen. As a result, PAH's in drinking water must be monitored.

US EPA method 550.1 is a method used by regulatory agencies to determine PAH's in drinking water by liquid-solid extraction (LSE) and high performance liquid chromatography (HPLC).<sup>1</sup> The method outlines the steps needed to perform extractions of PAHs from drinking water sources and finished drinking water; as well as the quality control (QC) program requirements to ensure the method is operated under control.

Automation of the extraction process allows more hands-off time for the operator and ensures better precision between operators. US EPA drinking water methods have embraced the use of SPE, yielding smaller extracts that require less evaporation time and require less solvent recovery. Use of Horizon Technology's automated sample preparation product line allows users to achieve the precise and consistent data needed to comply with EPA method 550.1 QC specifications as stated in method, section 10.3.2., for Initial Demonstration of Capability (IDC); while at the same time, streamlining laboratory practices to increase lab productivity.

#### Experimental

The instrumentation used in the method included the SPE-DEX<sup>®</sup> 5000 disk extraction system (Horizon Technology). This is a three-position extraction system that automates disk conditioning, sample loading and elution of the disk. Atlantic<sup>®</sup> C18 disks (Horizon Technology) were used for sorption of the analytes. The eluent was dried of residual water using DryDisk<sup>®</sup> membrane drying with the DryVap<sup>®</sup> In-line drying and evaporations system.

Analysis of the prepared extracts was accomplished with the Prominence i-Series HPLC (Shimadzu Scientific Instruments)



## Method Summary

1. Obtain 1-liter samples of DI water.
2. Acidify each sample with hydrochloric acid until pH <2.
3. Add surrogate compounds to the samples at a concentration of 5 µg/L.
4. Extract the samples using the SPE-DEX 5000 with 47-mm Atlantic C18 SPE Disks using the method shown in Table 1.
5. Dry and concentrate the extracts using the DryVap with a DryDisk to 0.5-mL final volume. See Table 2 for DryVap conditions.
6. Analyze the samples using conditions used in HPLC Table 3.

Table 1. SPE-DEX 5000 Extraction Method for 47-mm Atlantic C18 Disks

Step	Solvent	Solvent Vol. (mL)	Purge Time (s)	Pump Rate (s)	Sat. Time (s)	Soak Time (s)	Drain Time (s)	
1. Condition SPE Disk	MeCl	15	60	2	1	60	60	
2. Condition SPE Disk	MeOH	11	60	2	1	60	2	
3. Condition SPE Disk	Reagent Water	9	30	2	1	5	5	
4. Condition SPE Disk	Reagent Water	9	60	2	1	30	0	
Step		Sample Flow Rate (#)	Done Loading Sample Delay (s)					
5. Load Sample		2	45					
Step		Dry Time (s)	Pump Rate (#)	N2 Blanket				
6. Air Dry Disk Timer		60	6	Off				
Step	Solvent	Solvent Vol. (mL)	Purge Time (s)	Pump Rate (#)	N2 Blanket	Sat. Time (s)	Soak Time (s)	Elute Time (s)
7. Elute Sample Container	Acetonitrile	7	60	2	Off	1	60	45
8. Elute Sample Container	MeCl	7	15	2	Off	1	60	45
9. Elute Sample Container	MeCl	7	15	2	Off	1	60	45
10. Elute Sample Container	MeCl	7	15	2	Off	1	60	45
11. Elute Sample Container	MeCl	7	15	6	Off	2	60	60

Table 2. DryVap Conditions


Parameter	Setting
Dry Volume	20
Heat Power	5
Auto Rinse Mode	Off
Heat Timer	Off
Vacuum Pump	-8 psi
Nitrogen	20 psi

Table 3. HPLC Conditions

Parameter	Setting
HPLC System	Shimadzu i-Prominence Series
Column	Pinnacle® II PAH, 4 µm (250 mm x 4.6 mm) (Restek)
Detection	uV 220, 254 nm
Mobile Phase	60% ACN: 40% H <sub>2</sub> O
Elution	Gradient
Flow Rate	1.5 mL/min
Column Temperature	40 °C
Run Time	40 minutes

## Results and Discussion

Figure 1 shows a chromatogram of the analytes of interest in a standard mix. Resolution is good under the conditions specified.



ID #	Compound
1	Naphthalene
2	Acenaphthylene
3	Acenaphthene
4	Fluorene
5	Phenanthrene
6	Anthracene
7	Fluoranthene
8	Pyrene
9	Benz(a)anthracene
10	Benzo(k)fluoranthene
11	Benzo(b)fluoranthene
12	Chrysene
13	Benzo(a)pyrene
14	Dibenz(a,h)anthracene
15	Benzo(ghi)perylene
16	Indeno(1,2,3-cd)pyrene

Figure 1. Chromatogram of the 16 compounds of interest, identified in the key table.

Nine samples were extracted on the SPE-DEX 5000, spiked at a concentration of 5.0 µg/L. Three samples were run simultaneously, with each run taking 35-40 minutes to complete. Extracts were then dried and concentrated, taking approximately 1 hour to 1 hour 15 minutes to achieve a 0.5-mL final end volume. Final extracts were prepared and ready for HPLC analysis in under 2 hours.

The results for nine replicate samples are shown in Table 4. QC requirements set forth in EPA Method 550.1 Section 10.3, Initial Demonstration of Capability (IDC) state that the recovery value for at least three out of four consecutively analyzed samples must fall in the range of  $R \pm 30\%$  (or within  $R \pm 3$  Sr if broader) using the values of R (mean recovery) and Sr (standard deviation of the %) for reagent water in Table 2 of the method. Table 2 of EPA Method 550.1 shows compound recoveries ranging from 70.5 % - 93.6 %. While analyzing 9 samples using the SPE-DEX 5000, Atlantic C18 Disk and the DryVap, the mean compound recovery (R) was found to range from 80.2 % - 95.7 % recovery throughout each run. The relative standard deviation (Sr) of the mean recovery throughout these runs was calculated to be 5.1-6.3, meeting the method IDC criteria of  $\pm 30\%$ .

Table 4. Replicate Spiked Samples

Compound	%Rec.	% Rec.	S <sub>r</sub> (Std Dev)							
Naphthalene	71.9	77.0	82.3	73.0	81.6	91.7	80.1	84.7	88.3	6.2
Acenaphthylene	78.7	83.5	88.0	75.7	84.1	93.1	87.2	92.4	95.8	6.3
Acenaphthene	76.7	81.8	83.8	75.9	84.7	94.3	83.7	86.2	92.2	5.8
Fluorene	78.3	83.0	85.8	76.7	85.4	94.1	87.2	87.8	96.2	6.1
Phenanthrene	82.3	88.7	90.4	80.0	89.1	98.3	88.8	89.7	97.9	5.7
Anthracene	80.8	86.5	88.2	87.6	95.9	91.8	73.3	87.2	96.1	6.8
Fluoranthene	81.9	88.1	88.2	81.2	87.7	97.2	87.7	89.3	98.7	5.5
Pyrene	82.8	89.2	89.4	81.3	88.6	97.3	87.6	89.3	97.2	5.1
Benz(a)anthracene	82.1	88.3	88.0	81.6	87.0	96.9	85.8	89.0	99.3	5.6
Benzo(k)fluoranthene	82.1	88.3	87.9	82.1	86.9	96.7	85.0	89.2	99.7	5.7
Benzo(b)fluoranthene	82.5	88.8	88.9	82.9	88.0	97.7	86.4	89.8	99.9	5.6
Chrysene	82.0	88.4	87.6	82.2	86.6	95.7	83.7	87.8	99.2	5.5
Benzo(a)pyrene	81.0	87.5	87.8	81.7	86.7	94.4	80.6	86.7	97.8	5.5
Dibenz(a,h)anthracene	78.5	84.2	83.6	78.4	84.0	86.7	80.7	89.8	97.9	5.8
Benzo(ghi)perylene	80.6	87.3	87.2	81.5	85.9	93.9	82.9	88.2	99.0	5.6
Indeno(1,2,3-cd)pyrene	81.1	87.9	87.5	82.2	86.7	95.4	83.4	87.7	99.5	5.7
R (Mean Recovery)	80.2	86.2	87.2	80.2	86.8	95.7	84.0	88.4	97.2	

## Conclusion

Automated SPE for preparation of PAHs in drinking water gave excellent results, meeting the requirements set forth in US EPA method 550.1. In addition, the time for preparation was efficient and samples were prepared and ready for HPLC analysis in less than 2 hours. This is a powerful combination for laboratory efficiency.

## References:

1. US EPA Method 550.1, US EPA <https://www.o2si.com/docs/epa-method-550.1.pdf> (1990).

### EUROPE

Main Office: +46 18 565900  
Toll Free: +800 18 565710  
Fax: +46 18 591922  
Order Tel: +46 18 565710  
Order Fax: +46 18 565705  
order@biotage.com  
Support Tel: +46 18 56 59 11  
Support Fax: +46 18 56 57 11  
eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
Support Tel: +1 800 446 4752  
Outside US: +1 704 654 4900  
us-1-pointsupport@biotage.com

### HORIZON TECHNOLOGY

Tel: +1 603 893 3663  
Toll Free: +1 800 997 2997  
Fax: +1 603 893 4994  
sales@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
Fax: +86 21 68162829  
cn\_order@biotage.com  
cn-1-pointsupport@biotage.com

### KOREA

Tel: +82 31 706 8500  
Fax: +82 31 706 8510  
korea\_info@biotage.com  
kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 22 4005 3712  
india@biotage.com

Distributors in other regions are listed on [www.biotage.com](http://www.biotage.com)

### Part Number: AN119-HOR

© 2010 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.



a part of

