

Michael Ebitson, Horizon Technology, Inc.

Key Words

Endocrine Disrupting Hormones, drinking water, SPE, LC/MS/MS

APPLICATION NOTE



Introduction

Hormones in drinking water are a growing concern, since the presence of these compounds can potentially be linked to serious health hazards such as human developmental and reproductive _____

side-effects. Hormones are introduced into our environment in several ways. Natural and synthetic hormones that are given to livestock, pass though the animals, runoff into surface water, and leach into underground water supplies. Humans also produce and excrete natural hormone waste every day. An example is the natural and synthetic estrogens and androgens that are used in large quantities for medical purposes, such as oral contraceptives and hormone replacement. Human waste water flows through to sewage treatment plants, but these plants are not designed to handle the removal of hormones, so the hormones filter though the treatment process, and can ultimately end up in drinking water.

The Atlantic[®] HLB (Hydrophilic-Lipophilic Balanced) solid-phase extraction disk is able to retain the listed organic semi-volatile compounds through its unique ratio of hydrophilic N-vinylpyrrolidone and lipophilic divinylbenezene sorbent. The Atlantic HLB sorbent has a number of key features that enable it to perform the extraction, including stability at pH extremes and in a wide range of solvents, extraordinary retention of polar compounds, and a hydrophobic retention capacity three times greater than that of traditional silica-based SPE sorbents.

The SPE-DEX[®] 4790 provides automatic extraction of liquid samples by solid phase extraction methods, and can handle *Figure 1: Chemical Structures of Hormones used in this analysis.* samples ranging from 20 mL to 4 L. The Envision[®] Platform

Controller provides a user-friendly, web-based controller capable of interacting with up to eight extractors via one PC. The DryVap[®] Concentrator System provides automatic sample drying with patented DryDisk[®] separation membrane technology, automatically concentrating each dried extract by applying heat, vacuum, and sparge flow for up to six samples simultaneously. The Reclaimer[™] Solvent Recovery System (SRS) (now replaced by the SolventTrap _{SVOC}) is designed to handle the volume of solvent vapor generated from the DryVap Concentrator System while in operation.



Ø Biotage

Instrumentation

- Horizon Technology SPE-DEX 4790 Automated Extraction System
- Envision Controller
- DryVap Concentrator
- DryDisk Separation Membranes
- Reclaimer (Solvent Recovery System)
- Atlantic HLB SPE Disks (47mm)
- Agilent 1200SL HPLC
- Sciex API 4000 MS/MS

Method Summary

- 1. A 500 mL aliquot of sample is used.
- 2. Samples are preserved with 25 mg of ascorbic acid.
- 3. Agitate the bottles by hand until preservative is dissolved.
- 4. Verify the pH is between 5-7.
- 5. Spike the compound standard into samples (samples spiked at 50 ppt).
- 6. Start extraction method shown in Table 2 and programmed into the Envision Controller.
- 7. Collect the extract at high vacuum at -25 in. Hg.
- 8. Extract volume should be approximately 15 mL.
- 9. Add the extract to the DryDisk holder and start the concentration process on the DryVap using the conditions in Table 1.
- 10. Concentrate the extract to 0.5 mL.
- 11. Add 2 mL of MTBE into the concentrator tube and mix for solvent exchange.
- 12. Concentrate to 0.5 mL
- 13. Transfer the extract to a GC vial.
- 14. Analyze by LC/MS/MS.



Shown above: SPE-DEX 4790 Automated Extraction System, Envision Controller and DryVap Drying and Concentration System

Table 1: Conditions used for the DryVap Concentrator

| Parameter | Setting |
|------------|---------|
| Dry Volume | 20 |
| Heat Power | 5 |
| Auto Rinse | Off |
| Heat Timer | Off |

Table 2: Hormone Method Programmed into the Envision Controller

| Step | Solvent | Soak Time | Dry Time | |
|-------------------|---------------|-----------|----------|--|
| Prewet 1 | MTBE | 1:00 min | 30 sec | |
| Prewet 2 | Methanol | 1:00 min | 2 sec | |
| Prewet 3 | Reagent Water | 1:00 min | 0 min | |
| Prewet 4 | Reagent Water | 1:00 min | 0 min | |
| Process Sample | | | | |
| Air Dry 10:00 min | | | | |
| Rinse Step 1 | MTBE | 1:30 min | 30 sec | |
| Rinse Step 2 | MTBE | 1:30 min | 30 sec | |
| Rinse Step 3 | MTBE | 1:30 min | 1:30 min | |
| Rinse Step 4 | MTBE | 1:30 min | 2:00 min | |



Results

The results for the Atlantic HLB disk are listed in Table 3. The sample was spiked at 50 ppt and resulted in a concentration factor of 1000. The final extract concentration was 50 ppb. The table shows the compound name, amount recovered, and percent recovery. The Atlantic HLB SPE disk had a recovery range between 88-91% for the hormones. Table 4 shows the chromatograms of the selected hormones and their calibration curves.

Table 3: Recoveries of Selected Hormones

| Compound | Final Concentration (ppb) | Recovery (%) |
|-------------------|------------------------------|-----------------|
| Alpha-Estradiol | 45.4 | 91 |
| Ethynyl Estradiol | 44.2 | 88 |
| Progesterone | 44.3 | 89 |
| Esterone | 44.5 | 89 |
| Testosterone | 44.8 | 90 |

Table 4: Chromatograms and Calibration Curves for Selected Hormones





Analyte name: EthnylEstradiol (279.1 / 133.1) in blue





Analyte name: Progesterone (315.1 / 109.0) in blue



Analyte name: Progesterone (315.1 / 109.0) Calibration Equation: y = -0.09310 x^2 + 5708.08805 x + 10034.34686 (r = 0.99994)

Concentration













Conclusion

The resulting data from this study demonstrates that the Horizon Technology fully-automated extraction, drying and concentration systems used with the Atlantic HLB SPE disk are capable of the detection of selected hormones, resulting in data that is both accurate and precise. Extraction times were typically 30 minutes and drying and concentrating times were approximately 20 minutes.

The Horizon Technology SPE-DEX 4790 Automated Extractor System coupled with the Envision Platform, DryVap Concentrator System and the Reclaimer Solvent Recovery System reduces analyst labor, solvent usage, turn-around-time, and greatly improves accuracy and precision.

Acknowledgement

Horizon would like to thank Ali Haghani, from Mongomery Watson Labs (now Eurofins Eaton Analytical Laboratory) for his time and assistance with this study.

| www.horizontechinc.com | AN0011509_02 |
|--|---|
| 16 Northwestern Drive, Salem, NH 03079 USA | Tel: (603) 893-3663 Email: Support-Service@horizontechinc.com |

EUROPE

Main Office: +46 18 565900 Toll Free: +800 18 565710 Fax: +46 18 591922 Order Tel: +46 18 565710 Order Fax: +46 18 565705 order@biotage.com Support Tel: +46 18 56 59 11 Support Fax: +46 18 56 57 11 eu-1-pointsupport@biotage.com NORTH & LATIN AMERICA Main Office: +1 704 654 4900 Toll Free: +1 800 446 4752 Fax: +1 704 654 4910 Order Fa: +1 704 654 4900 Order Fax: +1 434 296 8217 ordermailbox@biotage.com Support Tel: +1 800 446 4752 Outside US: +1 704 654 4900 us-1-pointsupport@biotage.com HORIZON TECHNOLOGY Tel: +1 603 893 3663 Toll Free: +1 800 997 2997 Fax: +1 603 893 4994 sales@biotage.com

JAPAN Tel: +81 3 5627 3123 Fax: +81 3 5627 3121 jp_order@biotage.com jo-1-pointsupport@biotage.com CHINA Tel: +86 21 68162810 Fax: +86 21 68162829 cn_order@biotage.com cn-1-pointsupport@biotage.com

KOREA Tel: + 82 31 706 8500 Fax:+ 82 31 706 8510 korea_info@biotage.com kr-1-pointsupport@biotage.com INDIA Tel: +91 22 4005 3712 india@biotage.com

Distributors in other regions are listed on www.biotage.com



Part Number: AN001-HOR

© 2015 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE: A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.