

PN7140

FFF Solvent Organizer Module



PN7140 Solvent Organizer

Features

Operation Principle:

Holder for bottles and possibility of direct connection of up to 8 eluent lines via the integrated PEEK connector block.

Bottle Volumes:

For two 2 L, four 1 L, six 0.5 L or eight 0.25 L eluent/waste bottles.

Special versions may be provided by Postnova as well. Please inquire for possible options, which will depend on the requested customer application.

The PN7140 will be delivered including 2 x 1 liter bottles with bottle caps and 1/8" OD tubing.

Ordering Information

S-REA-7140-001 2-Channel Version

Extendable up to 8 channels



Specifications

- Working Principle:
By the use of a special collection tray
- Flow Rate:
Max. 10 mL/min for each channel
- Channels:
Standard: 2 channels
Extendable up to 8 channels
- Wetted Parts:
PEEK, ETFE, PFA, UHMW-PE
- Dimensions:
Width x Height x Length
270 x 90 x 430 mm
- Weight:
4 kg

Contact

- Postnova Analytics GmbH
86899 Landsberg, GERMANY
T: +49 8191 985 688 0
- Postnova Analytics UK Ltd.
Malvern, Worcestershire, WR14 3SZ, UK
T: +44 1684 585167
- Postnova Analytics Inc.
Salt Lake City, UT 84102, USA
T: +1 801 521 2004
- Postnova North Europe
01630 Vantaa, FINLAND
T: +358 9 8545 510

info@postnova.com
www.postnova.com