

ISOLUTE[®] Bulk Sorbents

ISOLUTE®

Bulk Sorbents

The same high quality sorbents used to manufacture ISOLUTE® SPE columns and 96-well plates are also available in bulk. The tightly controlled physical and chemical properties of these sorbents which provide reproducible performance in SPE procedures also ensure their suitability when used in applications that require loose material.

Non-Polar Sorbents

Part Number	Description	Qty.
9220-0025	ISOLUTE C18, Bulk	25 g
9220-0100	ISOLUTE C18, Bulk	100 g
9370-0100	ISOLUTE C18 MSPD, Bulk	100 g
9221-0025	ISOLUTE C18(EC)*, Bulk	25 g
9221-0100	ISOLUTE C18(EC)*, Bulk	100 g
9221-1000	ISOLUTE C18(EC)*, Bulk	1 Kg
9371-0100	ISOLUTE C18(EC)* MSPD, Bulk	100 g
9371-1000	ISOLUTE C18(EC)* MSPD, Bulk	1 Kg
9240-1000	ISOLUTE MFC18, Bulk	1 Kg
9290-0100	ISOLUTE C8, Bulk	100 g
9291-0100	ISOLUTE C8(EC)*, Bulk	100 g
9320-0100	ISOLUTE C2, Bulk	100 g
9321-0100	ISOLUTE C2(EC)*, Bulk	100 g

*(EC) - endcapped, a chemical process to reduce the concentration of silica surface silanol groups that provide polar and weak cation exchange secondary interactions. For more details see the QuickStart Guide to SPE. Literature part number UI331.

Mixed-Mode Sorbents

Part Number	Description	Qty.
9905-0010	ISOLUTE HCX-3, Bulk	10 g
9904-0025	ISOLUTE Multimode, Bulk	25 g

Polar Sorbents

Part Number	Description	Qty.
9460-0100	ISOLUTE SI, Bulk	100 g
9470-0025	ISOLUTE NH2, Bulk	25 g
9470-0100	ISOLUTE NH2, Bulk	100 g
9480-0100	ISOLUTE PSA, Bulk	100 g
9712-0100	ISOLUTE FL, Bulk	100 g
9430-0025	ISOLUTE DIOL, Bulk	25 g
9430-0100	ISOLUTE DIOL, Bulk	100 g
9714-0100	ISOLUTE AL-N, Bulk	100 g

Ion Exchange Sorbents

Part Number	Description	Qty.
9480-0100	ISOLUTE PSA, Bulk	100 g
9500-0025	ISOLUTE SAX, Bulk	25 g
9500-0100	ISOLUTE SAX, Bulk	100 g
9503-0100	ISOLUTE PE-AX, Bulk	100 g
9520-0100	ISOLUTE CBA, Bulk	100 g
9530-0025	ISOLUTE SCX, Bulk	25 g
9530-0100	ISOLUTE SCX, Bulk	100 g



ISOLUTE HM-N, a modified form of diatomaceous earth, is also available in bulk. The material is used for applications requiring an inert support as part of the sample preparation process (e.g. SFE). For details of pre-packed columns containing ISOLUTE HM-N, see [page 40](#).

Inert Support (modified diatomaceous earth)

Part Number	Description	Qty.
9800-0025	ISOLUTE HM-N, Bulk	25 g
9800-0500	ISOLUTE HM-N, Bulk	500 g
9800-1000	ISOLUTE HM-N, Bulk	1 Kg
9800-5000	ISOLUTE HM-N, Bulk	5 Kg
9800-10000	ISOLUTE HM-N, Bulk	10 Kg