





CAPSULE SINKERS

Distek Capsule Sinkers come in a wide variety of shapes and sizes to fit all your Dissolution Applications. The chart below highlights our complete offering including; Image, Part Number, Description, and Dimensions.


	PART NUMBER	DESCRIPTION	LENGTH	DIAMETER
	0500-0467	Sinker, Basket with Clip on Lid, 316 SS	26.9 mm	15.7 mm
	0500-0468	Sinker, Basket with Clip on Lid, 8 Mesh, 316 SS	20.1 mm	10.4 mm
	0500-0473	Sinker, Helical with Epoxy Coating, 6 Turn Spiral, Magnetic	21.3 mm	9.4 mm
	0500-0474	Sinker, 6 Turn Spiral, 316 SS	21.3 mm	9.4 mm
	0500-0475	Sinker, XL, Teflon Coated	29.2 mm	11.8 mm
	0500-0476	Sinker, "S" Shaped, 316 SS	38.1 mm	25.4 mm
	0500-0477	Sinker, 5 Turn, 316 SS	21.3 mm	9.8 mm
	0500-0478	Sinker, 4 Turn, SS, Inert Coated	21.3 mm	9.8 mm
	0500-0479	Sinker, 5 Turn, SS, Inert Coated	27.9 mm	10.4 mm
	0500-0491	Sinker, Swiss-Style, 316 SS	15.5 mm	5 mm
	0500-0492	Sinker, Swiss-Style, 316 SS	18 mm	6 mm
	0500-0493	Sinker, Swiss-Style, 316 SS	19.3 mm	7 mm

CAPSULE SINKERS

	PART NUMBER	DESCRIPTION	LENGTH	DIAMETER
	0500-0494	Sinker, Swiss-Style, 316 SS	23 mm	8 mm
	0500-0495	Sinker, Swiss-Style, 316 SS	31 mm	11 mm
	0500-0496	Sinker, Japanese Pharmacopeia Basket, 316 SS	25 - 26 mm	12 mm $\pm$.2
	5000-0011	Sinker, 3 Prong Black Plastic	25.4 mm	5.1 mm
	5000-0012	Sinker, 3 Prong Metallic Core, Inert Plastic Coated, Magnetic	25.4 mm	5.1 mm

COMMON CAPSULE SIZES

CAPSULES SIZE - CLOSED LENGTH



	SIZE 000	SIZE 00	SIZE 0EL	SIZE 0	SIZE 1	SIZE 2	SIZE 3	SIZE 4
MILLIMETERS	26	23.4	23.4	21.6	19.4	17.6	15.7	14.3
TOLERANCE	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3	$\pm$ 0.3
INCHES	1.024	0.921	0.921	0.85	0.764	0.693	0.618	0.563
TOLERANCE	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012	$\pm$ 0.012