

## BIOROCESS CONTROL STATION



## THE NEW interFACE of **BIOREACTORS!**

Distek, Inc. has revolutionized benchtop scale bioreactors with the release of our BIOne bioprocess control station for both mammalian and microbial models. The sleek appearance and intuitive touchscreen interface are only the beginning. The BIOne controller can seamlessly transition between traditional glass reactors and our BIOne single-use system. Easy to use, the BIOne bioprocess control station is powerful enough to develop processes at the benchtop scale. Designed to fit your needs, the BIOne controller can be customized to conform to your lab design... rotameters or auto flow control, left or right side setup, optical pH or standard, smart sensors or standard sensors and much more.



BIONQ

FEATURES, ADVANTAGES & BENEFITS

Touchscreen Interface

Smart Media Filling

Email and Text Alerts

• Trend up to 8 Parameters

• Variable Speed Pumps

• Available with Single-Use or Glass **Bioreactors** 

• 2L, 5L & 10L Working Volumes (10L Available in Glass only)

Configurable for your Lab • Up to 5 Gas Control Modules • Rotameters or Auto Flow Control • Optical pH or Standard • Left or Right Side Setup





121 North Center Drive North Brunswick, NJ 08902 888.2.DISTEK distekinc.com bione@distekinc.com



Distek, Inc. is a leading manufacturer of laboratory testing instruments for the pharmaceutical and biotechnology industry for over 40 years, as well as an experienced provider of validation and qualification services. Distek engineers highly innovative, technically modern and user-friendly instruments that provide many innovative features to customers, worldwide. Distek's bioprocessing portfolio includes the BIOne single-use bioreactor for cell culture applications, as well as an innovative bioreactor control system for both mammalian and microbial models. Distek's core market includes brand name and biosimilar manufacturers, CRO's, CMO's, government agencies and universities. Distek holds several instrumentation and accessory patents and has been ISO certified since 2002. For technological advancements that offer reliable and consistent performance, the biotech industry trusts Distek.

